


Blue Heron pUC119



**Blue Heron Bio pUC
(3150 bp)**

(pUC119 derivative)

polylinker

amp resistant

bi-directional

You may download the [Blue Heron Bio pUC map](#) which shows the sequence of the adaptors, and/or the [Blue Heron Bio pUC sequence](#) file.

Blue Heron Bio may flank your gene with up to 50 bp unique sequence to facilitate Blue Heron Bio cloning. Your gene and any unique flanking sequence is designed by Blue Heron Bio to to be fused non-directionally into the Blue Heron Bio adaptors between the CTTTC and GAAAG "sticky ends". If a unique flanking sequence is used, a cloning information sheet will be posted in your File Area: Notes from Blue Heron folder and included with your packing slip when your gene is shipped. This unique flanking sequence will not contain common restriction sites.


M13/pUC Forward 23-Base Sequencing Primer: 5'-CCCAG TCACG ACGTT GTAAA ACG-3'

M13/pUC Reverse 23-Base Sequencing Primer: 5'-AGCGG ATAAC AATTT CACAC AGG-3'

BLUEHERON

BIOTECHNOLOGY

Blue Heron pUC



Sequence of polylinker flanking the insert

```
GAATTCATGGTCTCAACTTTC(***)Synthetic gene by GeneMaker™(***)GAAAGGAGATAGGATCCAAGCTT  
EcoRI NcoI BamHI HindIII
```

*** Your gene may be flanked by up to 18 bp to facilitate Blue Heron cloning. This will not contain any common restriction sites. If your gene is cloned with this sequence a cloning information sheet will be included with your gene shipment and electronically posted on your GeneMaker Secure Website File Area in the "Notes from Blue Heron" folder. We recommend you frame your gene with restriction sites to facilitate your future subcloning.

Blue Heron pUC119 DNA Sequence

```
1 agcgccaat acgcaaaccg cctctccccg cgcggtggcc gattcattaa tgcagctggc acgacaggtt tcccgactgg aaagcgggca gtgagcgcaa
101 cgcaattaat gtgagttagc tcaactcatta ggcaccccag gctttacact t tatgcttcc ggctcgtatg ttgtgtggaa ttgtgagcgg ataacaattt
201 cacacaggaa acagctatga ccatgattac gccaaagcttg catgcoctgca ggtcgaactc agaggatccc cgggtaccga gctcgaattc actggccgtc
301 gttttacaac gtcgtgactg ggaaaaccct ggcgttaccg aacttaactg ccttgcaaga catccccctt tcgccagctg gcgtaatagc gaagaggccc
401 gcaccgatcg ccttcccaa cagttgcgca gcctgaatgg cgaatggcgc ctgatgcggt attttctcct tacgcatctg tgcggtattt cacaccgat
501 acgtcaaagc aaccatagta cgcgcctgt agcggcgcat taagcggcgc ggggtgtggg gttacgcgca gcgtgaccgc tacacttgc agcgccttag
601 cgcccgtccc tttcgctttc ttcccttcc tctcgccac gttcgccggc tttcccgtc aagctctaaa tcgggggctc ctttagggg tccgatttag
701 tgctttacgg cacctcgacc ccaaaaaact tgat ttgggt gatggttcac gtagtgggc atcgcctga tagacggtt ttcgccctt gacgttgag
801 tccacgttct ttaatagtgg actcttgttc caaactggaa caaactcaa ccctatctcg ggctattctt ttgattata agggattttg ccgattcgg
901 cctattggtt aaaaaatgag ctgatttaac aaaaatthaa cgcgaatttt aacaaaatat taacgtttac aattttatgg tgcactctca gtacaactg
1001 ctctgatgcc gcatagttaa gccagccccg acaccggcca acaccgctg acgcgcctg acgggcttgt ctgctcccgg catccgctta cagacaagct
1101 gtgaccgtct cggggagctg catgtgtcag aggttttcac cgtcatcacc gaaacgcgcg agacgaaagg gcctcgtgat acgcctattt ttaggtta
1201 atgtcatgat aataatggtt tcttagacgt cagggtggc ttttcgggga aatgtgcgcg gaaccctat ttgtttat ttctaaatac attcaaat
1301 gtatccgctc atgagacaat aaccctgata aatgcttcaa taatattgaa aaaggaagag tatgagtatt caacatttcc gtgtcgccct tattccctt
1401 tttgcccgat tttgccttcc tgtttttgct caccagaaa cgctgggtgaa agtaaaagat gctgaagatc agttgggtgc acgagtgggt tacatcgaa
1501 tggatctcaa cagcggtaag atccttgaga gttttcgccc cgaagaacgt tttccaatga tgagcacttt taaagtctg ctatgtggcg cggattatc
1601 ccgtattgac gccgggcaag agcaactcgg tcgcccata cactattctc agaatgactt ggttgagtac tcaccagtca cagaaaagca tcttacgat
1701 ggcatgacag taagagaatt atgcagtgt gccataacca tgagtataa cactgcggcc aacttacttc tgacaacgat cggaggaccg aaggagctaa
1801 ccgctttttt gcacaaatg ggggatcatg taactcgctc tgatcgttgg gaaccggagc tgaatgaagc cataccaaac gacgagcgtg acaccacgat
1901 gcctgtagca atggcaaaa cgttcgcaa actattaact ggcgaactac ttactctagc ttcccggcaa caat taatag actggatgga ggcggataaa
2001 gttgcaggac cacttctgcg ctccggccctt ccggctggct ggtttattgc tgataaatct ggagcgggtg agcgtgggtc tcgcggtatc attgcagcac
2101 tggggccaga tggtaaagccc tcccgtatcg tagttatcta cacgacgggg agtcaggcaa ctatggatga acgaaataga cagatcgctg agataggtgc
2201 ctactgatt aagcattggt aactgtcaga ccaagtttac tcatatatac tttagattga tttaaaactt catttttaat ttaaaaggat ctaggatgag
2301 atcctttttg ataactcat gaccaaaatc ccttaacgtg agt tttcgtt cactgagcg tcagaccccg tagaaaagat caaaggatct tcttgagatc
2401 ctttttttct gcgcgtaatc tgctgcttgc aaacaaaaaa accaccgcta ccagcgggtg tttgtttgcc ggatcaagag ctaccaactc tttttccgaa
2501 ggt aactggc ttcagcagag cgcagatacc aaat actgtc cttctagtgt agccgtagt aggccaccac ttcaagaact ctgtagcac gcctacatac
2601 ctgctctgct taatcctggt accagtggct gctgccagt gcgataagtc gtgtcttacc gggttggact caagacgata gttaccgat aaggcgcagc
2701 ggtcgggctg aacgggggggt tcgtgcacac agcccagctt ggagcgaacg acctacaccg aactgagata ccta cagcgt gagctatgag aaagcgcac
2801 gct tcccga gggagaaagg cggacaggta tccggtaagc ggcagggtcg gaacaggaga gcgca cgagg gagcttccag ggggaaacgc ctggtatctt
2901 tatagctctg tcgggtttcg ccacctgta cttgagcgtc gat ttttgtg atgctcgtca ggggggcgga gcctatggaa aaacgccagc aacgcggcct
3001 ttttaacggtt cctggccttt tgctggcctt ttgctcacaat gttctttcct cgttatccc ctgatctgt ggataaccgt attaccgctt ttgagtgagc
3101 tgataccgct cgcccagacc gaacgaccga gcgcagcgag tcagtgagcg aggaagcgga ag
```

Blue Heron pUC119 Restriction Map

RESTRICTION MAP: ALL ENZYMES

```

=====
Enzyme      Cuts    [ ----- Base 5' to Cleave Sites ----- ]
=====
AatII       1      1230
AccI        1       253
AclI        3       963  1547  1920
AcyI        3       447  1227  1609
AflII       0
AflIII      1      3037
AgeI        0
AhaIII      3      1571  2263  2282
AluI        18       55  119  214  236  281  377  663  920
           1098 1117 1796 1859 1959 2480 2737 2783
           2873 3099
AlwNI       1      2628
ApaBI       0
ApaI        0
ApaLI       3       980  1477  2723
ApoI        3       285  933  944
AscI        0
AsuI        7       396  756  1169  1785  2007  2024  2103
  
```

AsuII	0										
AvaI	2	269	857								
AvaII	2	1785	2007								
AvrII	0										
BalI	0										
BamHI	1	264									
BclI	0										
BetI	3	1853	2684	2831							
BglI	2	437	2031								
BglII	0										
BsaAI	1	750									
BsaBI	0										
BseMII	2	1001	1662								
BsePI	0										
BsiYI	9	34	531	857	858	1111	2559	2838	3004	3022	
Bsp1407I	0										
BspHI	3	1204	1309	2317							
BspLU11I	1	3037									
BspMII	0										
BstEII	0										
BstXI	0										
Cac8I	17	57	86	238	242	246	375	379	399		
		590	604	647	1023	2036	2427	2987	3024	3110	
CauII	7	270	271	1077	1112	1613	1964	2660			
Cfr10I	2	645	2064								
CfrI	3	36	293	1756							

HaeI	3	2563	3015	3026					
HaeII	5	5	450	595	603	2797			
HaeIII	13	38	295	397	758	900	1171	1758	2025
		2105	2563	2997	3015	3026			
HgiAI	5	283	984	1481	1566	2727			
HgiCI	5	131	273	446	709	2196			
HgiJII	2	283	679						
HhaI	23	4	32	97	428	449	524	537	546
		568	594	602	1055	1158	1258	1590	1927
		2020	2413	2522	2696	2796	2863	3133	
HincII	1	254							
HindIII	1	234							
HinfI	8	41	255	798	820	2150	2667	3063	3138
HpaI	0								
HpaII	14	159	270	646	1077	1111	1612	1854	1964
		2031	2065	2469	2659	2685	2832		
KpnI	1	277							
MaeI	5	259	597	1956	2291	2544			
MaeII	11	309	501	639	749	792	804	963	1227
		1547	1920	2336					
MaeIII	13	313	333	560	572	1100	1488	1676	1829
		1887	2218	2501	2617	2680			
MboI	15	264	405	1466	1502	1519	1777	1823	1841
		2182	2287	2299	2377	2385	2396	2471	
McrI	5	408	1631	1780	2703	3127			
MfeI	0								
MluI	0								
MseI	20	47	106	344	540	811	909	926	937
		949	960	1017	1198	1570	1935	1974	2209
		2262	2276	2281	2333				

MslI	3	1378	1737	1896						
MstI	2	427	1926							
MwoI	17	10	94	138	237	407	437	446	453	
		541	543	585	612	642	2031	2419	2991	3105
NaeI	1	647								
NarI	1	447								
NcoI	0									
NdeI	0									
NheI	0									
NlaIII	11	225	244	1124	1208	1313	1706	1742	1820	
		1830	2321	3041						
NlaIV	14	133	266	275	448	678	690	711	1262	
		1852	2063	2104	2198	2970	3009			
NotI	0									
NruI	0									
NspBII	6	55	377	1047	1513	2454	2699			
NspI	3	244	1124	3041						
PacI	0									
PflMI	0									
PmaCI	0									
PmeI	0									
PpuMI	0									
PshAI	0									
PstI	1	250								
PvuI	2	408	1780							
PvuII	2	55	377							

SspI	2	958	1344							
StuI	0									
StyI	0									
SwaI	0									
TaqI	5	253	283	715	1495	2939				
TatI	2	990	1666							
TfiI	2	41	3063							
TseI	13	52	355	428	568	1117	1728	2094	2422	
		2628	2631	2696	3115	3133				
Tsp45I	5	313	572	1100	1676	1887				
Tsp4CI	9	422	776	1106	1141	1709	2224	2537	3007	3078
TspEI	10	103	178	195	285	933	944	970	1716	
		1971	2277							
TspRI	10	95	296	1730	1757	2104	2209	2358	2629	
		2642	3148							
Tth111I	0									
VspI	3	47	106	1974						
XbaI	1	258								
XcmI	0									
XhoI	0									
XhoII	7	264	1502	1519	2287	2299	2385	2396		
XmaIII	0									
XmnI	1	1549								

=====